

Screen description

Memory function



Recollect memorized data

Show energy level

Emmission of beam per second



Stop Cooling



Start Cooling

Hair Removal Quick-Reference Guide

AREA	FLUENCE	PASS
armpit/ Bikini line	FHR(20~30J) HR(30~40J,3Hz)	6 ~ 8 1 ~ 3
Arm & leg	FHR(30~40J) HR(40~50J, 3Hz)	6 ~ 8 1 ~ 3
Philtrum & chin	FHR(20~30J) HR(30~40J,3Hz)	6 ~ 8 1 ~ 3
forehead	FHR(20 ~ 30J) HR(30~40J,3Hz)	6 ~ 8 1 ~ 3
REJUVENATION	25J~35J 5Hz	A Part 15Minute

This parameter is a guideline to help users' treatment, and it may vary slightly according to patient's condition and surroundings.

How to work the Machine



Connect handpiece connector by pushing both round buttons of it. After connecting, make sure it is connected well by pulling it slightly.



Turn on the switch key



Make sure LCD screen shows up normally



Decide treating area and Select mode and adjust Jouls and Herz referring to parameter.



Pick up the handpiece and put on the strap band on the wrist .



Start treatment.

The manufacturer and distributor don't take responsibilities when damage of machine happened while not wearing strap band.

Error description

Aroma was designed to maximize safety for both patient and user. The below screens show up on LCD, refer to the measures.



Interlock check error. Make sure that interlock switch is fastened in the back below




Handpiece connection error
Make sure that hand piece is connected well in the back below .



It happens when the cooling water is not circulated well.
Turn off the power and get help from distributor.


CAUTIONS

The user should operate the machine after being fully aware of manual.



After treatment, be sure to wipe siphon tip with a wet towel.

- ▶ (If it is used with fur or debris attached on the tip, the patient can be burned and the tip can also be damaged.)



How to put in distilled water

Loosen screws and take off the panel

- ▶ Fasten funnel to inlet (any one of two inlets), put in distilled water to 90% of water container. (You can see the water filled through transparent column of water container.)

*After turning on, wait 10 minutes to warm up.

*After using 30 minutes, 10 minutes' pause is recommended.

*When the user doesn't use the machine for a long time, turn off the power.

*The user must put in distilled (or deionized) water, or module can be damaged.

*If the user disassembles and reassembles the machine without consent of the manufacturer and doesn't use distilled water, he (the user) won't get warranties.